Neuromuskulärer Stimulator
Neuromuscular Stimulator

Rehalito

Medizinprodukt / Medical Device

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Content

Please read the Operators Manual carefully
before using the device and pay attention
to all references !

	Page
I. Musclemulator Rehalito	
Delivery & Accessories / Medical Device Reporting	24
Indication / Contraindication	25
Warnings	26
Advice on Application of Electrodes	27
Technical Data	28
Explanation of Symbols / Maintenance	29
Warranty / Safety Adjustments	30
Real Time Clock Setup	31
II. Display Elements of the Rehalito	
Display and Operating Elements	32
Description of the Device	33
III. Application of the Rehalito	
Preset Programs	34-35
Setting the Custom Program	36-37
IV. Programs 01-05	38
Abbreviations	39
V. Lock-Button and Statistics	40-41
Trouble Shooting Guide	42
Placement of Electrodes	43-57
Waste Disposal / Recycling of Batteries / Literature / Copyright	58
Anatomy	59
Electromagnetic Compatibility	63-65

Delivery & Accessories

Included in the delivery:

REF: 010E-201: 1 device **Rehalito**, 2 lead wires, 4 self adhesive electrodes (40 x 80 mm), one 9V battery and 1 operators manual

With regular average use, the expected service life of **Rehalito** is 10 years. Service life of accessories is one year.

Accessories:

REF

011E-006: 1 set self adhesive electrodes 40x40 (4 pieces)

011E-007: 1 set self adhesive electrodes 40x80 (4 pieces)

011E-051: 1 set self adhesive electrodes 32mm Ø (4 pieces)

Spare parts:

REF:

012E-012: 2 lead wires

052F-804: 1 battery 9V

012E-209: 1 operators manual **Rehalito**

Rehalito shall be used with original accessories only.

The battery life is approximately 12 hours or more at 20 mA per channel.

The service life (expiry date) of electrodes and electrode jelly is printed on each package.

Medical Device Reporting:

If there is a serious deterioration in the state of your health when using the medical device, report this to your distributor or to the manufacturer and to the competent authorities of your country.

Indications

Rehalito was designed for the following indications and is easy to use without any assistance:

Muscle Stimulation:

Improvement of blood flow (P 3)

Improvement of muscle spasticity (P 3)

Muscle Strengthening

- Avoidance of muscular imbalances (P 1, 2)
- Activation of muscle fibers (P 1, 2)
- Muscle Fatigue (P 1)
- Strengthening of atrophied muscles (P 1, 2)
- Muscle Recovery (P 1, 2)

Conservation and extension of mobility (P 1, 2)

Nerve Stimulation:

Increase of release of bodily endorphins

General Pain Relief and Pain Therapy including but not limited to:

- Arthritis
- Back Pain including but not limited to:
 - Chronic Back Pain
 - Cervical Spine Pain
 - Thoracic Spine Pain
 - Lumbar Spine Pain
- Circulation Problems
- Neuralgia
- Myalgia
- Phantom Pain
- Fracture Pain
- Pelvic Floor Pain

For all these indications you can use program 4.

Contraindications

Rehalito may not be used in the following cases:

By patients with electronic implants like pumps or a pace maker because functional disturbances can occur.

By patients with myocardial damage (damage to the cardiac muscle) or arrhythmia (inconsistent heart beat).

In presence of flammable anesthetic mixtures with air or with oxygen or nitrous oxide.

When operating machinery with a potential of danger.

While driving a vehicle.

By patients with epilepsy.

If severe local inflammation or thrombosis is present.

By patients with mechanical ileus.

By patients with febrile illnesses or infectious diseases.

Do not use with a surgical high frequency apparatus because burn marks could occur underneath the electrodes.

When using metal-containing contraceptive spirals.

Warnings

Please pay unconditional attention to the following basic rules:

- Do not stimulate transcerebral (via the head).
- Do not stimulate via the Carotis-Sinus-Nerve (carotid artery).
- Do not stimulate in the area of the Sinus knot (heart).
- Do not stimulate on the Bulbous (eye).
- Do not stimulate in the Larynx-Pharynx-Area (throat/mouth).
- Do not stimulate in places with skin defects or broken skin.
- Do not position the electrodes directly on the spine.

Take special care with the use of the electric stimulation:

- if after an operation the healing process can be interrupted through a muscle contraction.
- after acute bleeding or bone fractures with a tendency to bleed.
- during menstruation or pregnancy.
- in case of desensitized skin with nerve damage.
- if stimulation is located close to osteosynthesis material (metal in the body).
- if patients cannot use the device properly, due to mental or physical limitations.

Please note also the following instructions:

Rehalito is type BF equipment and may not be immersed in liquids.

If device operation occurs close to short wave, microwave, cellular or other wireless equipment, instability of the values of the stimulator could occur causing a temporary reduction in current. The application of the electrodes nearby the thorax can increase the risk of cardiac fibrillation. Patients with metal implants and/or electronic may only stimulate with permission of an experienced doctor.

Stimulation in the areas of tumors is allowed only by instruction of the specialized physician.

The **Rehalito** may not be operated if the battery compartment is not closed.

Before opening the battery compartment, make sure that there are no lead wires in the jacks at the head of the unit to eliminate any possibility of current flow.

The operator should not touch the battery contacts & patient simultaneously.

Never connect the device directly to a battery charger or to any other mains powered equipment.

If using rechargeable batteries be sure to use a CE approved battery charger.

If you have swallowed any part of the device or any of its accessories, please contact a doctor immediately.

TENS and muscle stimulation should be used during menstruation or by pregnant women as well as by children under the guidance of a doctor or therapist only.

Remove batteries and battery cover, drain water from device and allow (Rehalito) to dry normally for 48 hours if accidentally submerged in water.

Please, keep this device stored out of the reach of children. Do not leave the lead wires in reach of children, since strangulation is possible due to the length of the lead wires.

Advice on Application of Electrodes

Make sure that **Rehalito** is turned off (no indication visible on the LCD display).

Now connect the single poled ends of the lead wires with the electrodes.

The output plugs of the electrode lead wires will be plugged into the channels A or B of the **Rehalito**.

Make sure, that all plugs fit completely into the electrode sockets as well as the device sockets. Please don't use any oil, ointment or cream, since all these have an isolating effect.

Now place the electrodes according to the instructions of your doctor or physiotherapist.

In the case of the skin gel losing its adhesion, you can't achieve more than 6 mA - please replace the electrodes. The same effect, however, also occurs when using a damaged lead wire.

Therefore, if the normal treatment is not possible by using new electrodes, the lead wire may be faulty.

Electrodes which were damaged due to an improper handling may not be used.

Electrodes with a damaged lead wire isolation (e.g. visible copper wire) may not be used.

Electrodes with damaged, missing or stained aluminium coating or carbon layer, also with insufficient adhesive gel may not be used.

Electrodes may not be altered (e.g. by means of scissors or other tools).

Do not pull the electrodes by force on the electrode lead wire.

Apply the electrodes on clean and healthy skin only.

In the case of strong bodily hair growth, the hair should be shortened with scissors.

Neither stretch the skin nor the electrodes when applying the electrodes.

To remove the electrodes from the skin, please lift the electrodes at the rim towards the middle.

Effective current densities of more than 2 mA/cm² require an increased attention of the user with regard to the intensity adjustment.

If you use self-adhesive electrodes, we recommend the use of our self-adhesive electrodes with an area of 40x80 mm.

Electrodes can contain amounts of nickel. Allergy sufferers should check this before using.

Should skin irritation occur whilst using the **Rehalito**, please interrupt the application and inform your doctor - to exclude an allergic reaction.

All electrodes are made for use of one single person only.

Technical Data

Channels:	Two isolated independent channels
Waveform:	Asymmetrical, rectangular, bi-phasic with zero DC current
Intensity:	0 - 80 mA into 500 ohm load 80 mA maximum into 1000 Ohm load 70 mA maximum into 1500 Ohm load (Indication only, actual intensity varies depending on condition of skin electrodes)
Classification:	Internally powered, constant current, continuous operation, Type BF applied part, class IIa, IP20
Frequency:	1 - 150 Hz (Burst 35 - 150 Hz)
Pulse Width:	50 - 450 μ S (Burst 50 - 200 μ S)
Time Settings:	Variable 1 minute to 9 hours 59 minutes (custom program)
Max. Voltage:	180 V without load (70V with 1000 ohm load)
Battery;	9 Volt 6LR61 / 6AM6 recommended

Dimensions of the device: H 11,7 cm / W 6,6 cm / D 3,3 cm

Weight of the unit: 150 g (incl. battery)

Firmware version: 3.1 or later

Max. environmental conditions for transportation: -25 to +70 degree centigrade, 0-93% atmospheric humidity. Recommended environmental conditions for storage: -25 to +40 degree centigrade, 0-93% atmospheric humidity. The device shall not be stored in direct sunlight. Allow up to one hour for the device to return to normal operating temperature if stored in conditions below +5°C or above +40°C.

Environmental conditions for use: +5 to +40 degree centigrade, 0-93% atmospheric humidity, air pressure 700-1060 hPa. During use the surface temperature of the device can be 41°C or higher, but never exceeds 48°C.

The **Rehalito** is protected against debris larger than 12.5 mm in diameter from entering the cabinet, is not water resistant and rated IP20.

The patient is the intended operator for this unit. This device is quality controlled and fulfills all legal requirements for medical products. Therefore this device carries the CE-sign „CE 0123“. The serial number is on the rear of the unit.

Manufacturer: MTR+ Vertriebs GmbH,
Kamenzer Damm 78, 12249 Berlin, Germany
www.mtrplus.com

Explanation of Symbols

	Attention Warning		Manufacturer
	Typ BF applied parts		Date of manufacture 2022 - 04
	Serial-Number of the product		CE-mark: conforms to essential requirements of the European standards
	Keep dry		Not suited for persons with cardiac pacemaker
	Batch Code		Read the instructions
	Medical Device		Unique Device Identification
	Temperature limitation (transport)		Keep away from sunlight
	Temperature limitation (storage)		Humidity limitation

Maintenance

Rehalito has been designed to be maintenance free. Nevertheless, please adhere to the following advice:

If needed, please clean **Rehalito** with a mild detergent, water and a damp cloth. Do not expose **Rehalito** to extreme humidity and do not submerge under water. You may clean the lead wires with a damp cloth. You may also from time to time treat the lead wires with talcum powder in order to prevent cracks and to improve the lead wire life.

A clean cloth, moistened with an alcohol free disinfection liquid can be used to clean the device. Please follow all instructions of the manufacturer of the liquid.

After termination of the therapy, **Rehalito** can be re-used by other persons. Please clean it before use with a mild disinfectant cleaner.

When changing the battery, please make sure to insert the new battery within 30 minutes after removal of the empty battery to ensure the devices time and date is kept current. We recommend that you always have a spare battery with the device.

The **Rehalito** does not require calibration.

The **Rehalito** is not repairable and contains no user serveable parts.

No modification to the **Rehalito** is allowed.

Warranty

MTR+ Vertriebs GmbH provides a 24 month warranty, commencing from the date of our delivery (invoice date), that in case of failure, the **Rehalito** will be serviced or replaced free of charge. All returns must first be authorised by MTR+ Vertriebs GmbH in advance.


Excluded from this warranty are failures of the **Rehalito** which result from excessive use, careless or improper handling as well as damage by force.

The warranty loses its lawful meaning if the producer MTR+ Vertriebs GmbH has not given its consent in writing prior to repairworks or servicing being carried out. Without this written consent the manufacturer will not return the device.

The manufacturer reserves the right of technical changes, which may serve to increase treatment success, the functionality of the device or its accessories.

Safety Adjustments

1. The +buttons become inactive after one minute running time of a program to avoid an unintentional increase of intensity . This lock is switched off again by pressing the +button three-times repeatedly. An increasing of intensity is possible again for another minute.
2. If the electrical circuit should be interrupted during operation (> 6 mA), for example by losing an electrode or by a broken lead wire, the device stops stimulating immediately and the display shows the flashing symbol of an electrode. The effected channel will also flash.
3. If any phase comes to its end the intensity falls to 8 mA automatically and must be increased by the user to the desired intensity once again.

Real Time Clock Setup

Rehalito has a “Real Time Clock”, that shows year, month, day, hours and minutes and allows to store detailed statistics of treatment data.

To set time and date, please select the custom program 05 and press the **S**-button for at least 5 seconds. You’ll reach the configuration-mode. Now press **S** once again for 10 seconds.

With the +/- button of channel A (left side) you can step forward to the following values

- hour (indicated by: H)
- minute (indicated by: MI)
- day (indicated by: dA)
- month (indicated by: MO)
- year (indicated by: Y)

The respective value is flashing and can be changed with the +/- buttons of channel B (right side). Finally you save all your settings by pressing the PRG-button once again.

Before starting the therapy, please set the clock once if necessary.

If **Rehalito** has been without battery for more than 30 minutes, the clock may need to be reset.

Display and Operating Elements

„On/Off“: By pressing the button once, the device will either be switched on or off.

„+/- buttons“: These buttons adjust the intensity for channels A and B. The device has two independently controllable output channels A and B. When pressing the buttons the intensity will either be increased or decreased in steps of 1 mA. 60 seconds after the last change of intensity, a further increase may not be carried out due to an electronic barrier. Pressing the „+“ button three times releases the barrier and the intensity can be altered again.

A reduction of the intensity is possible at any time.

Whilst configuring program PC05 you may switch between the programmable parameters by pressing the left „+/-“ button (channel A), and through pressing the right „+/-“ button (channel B) the desired value/data for respective parameter can be set. This way, the application time, the stimulation-mode, the frequency and the pulse width can be selected.

„P“: This button selects preset programs 01–04 or custom program PC05.

Press the „P“-button once at any time during a running program to reduce the mA on both channels to zero and pause the program. Pressing „P“ a second time ends the program.

Holding „P“ for three seconds while statistics are displayed will exit the statistics screen without deleting them. **Rehalito** will return to the home screen of the selected program and the lock will be cancelled.

The „P“ button stores all adjustments made to the custom program mode and returns to the home screen.

„S“: This button is used to setup custom program 05. Press the „S“-button for at least 3 seconds in program PC05 to enter the custom program setting. Pressing „S“ a second time for 10 seconds will enter the real time clock (RTC) setup.



Press the „S“ button once at any time during a running program, to reduce the mA on both channels to zero and pause the program. Pressing „S“ a second time will resume the program and the current will increase over a period of one second back to the previously set intensity level.

Holding the „S“ button for 5 seconds while the statistics are displayed, will display „dEL“.

Holding „S“ a second time for 5 seconds will delete the statistics, set the custom program to default values and show the RTC setup screen. Turn off the unit if you do not want to delete the statistics.

The unit is intended to be supported by hand during normal use. The belt clip is for temporary relief only.

Description of the Device

Preset Programs (01 - 04)

Preparations:

Make sure that **Rehalito** is switched off. Cleanse the area of the skin which is to be treated with clean water. Following that, make sure that all plugs are completely connected to the respective electrode sockets as well as device sockets. Place the adhesive electrodes onto the area that has been cleansed, according to the advices in pages 23-25 of this manual.

Rehalito is operated by a 9 volt block battery. Open the battery compartment at the back by pushing down the lid and insert the battery by hand. Please ensure the correct polarity in accordance with the markings. **Rehalito** is now operational.

If the battery voltage is too low, the battery-symbol in the LCD-display will flash. Please replace the battery at this time. Empty batteries are special waste and therefore have to be disposed off at an appropriate collection point or recycled. Never dispose of batteries via domestic or household waste.


Adjustment:

1. Press the button „On/Off” once.

Rehalito is now switched on. The program used last, appears on the LCD-display. Program parameters will be displayed as well as the treatment mode and the total treatment time.

2. By pressing the program button „P” you can choose from the programs 01 – 04 (pre-programmed).
3. This is followed by pressing the buttons „+/-” below the channel used. The selected program runs automatically and the LCD displays the intensity as milliamps (mA).

In the upper row you can see the remaining program time and in the bottom line the intensity being applied to each channel. If the button „+” is pressed repeatedly the intensity of stimulation increases. If the button „-” is pressed repeatedly the intensity of stimulation is reduced. For the setting of the intensity it is vital that one feels the electric current but that discomfort never occurs.

4. If you wish to interrupt the program, press the button „S” once and the program will pause. The program will resume for the remaining program time when pressing this button once again. In this case, the intensity increases within one second to the previous used value.
5. By pressing the „P”-button you pause as well, you can continue with new intensity adjustment by pressing the +button. If you press „P” a second time the program is finished.
6. **Rehalito** is switched off when pressing the button „On/Off” or automatically one minute after the program is finished.

Setting the Custom Program (05)

1. Switch on **Rehalito** by pressing the „**On/Off**” button.
2. Choose the program PC05 with the „**P**” button.
3. Press the „**S**” button for at least three seconds. **Rehalito** changes to the programming mode and the time that indicates program duration, flashes to show that this function can be altered.

In the following please use the left „+“ button to store your choice and to change to the next step. The right „+/-“ buttons are respectively used in order to alter the chosen and flashing parameter.

4. Select with the right „+/-”buttons the desired time for the phase.
5. By pressing the left „+”button you will be able to store your choice and change to the next step. Mode now flashes on the LCD and you are able to select the mode that you require by pressing the right „+/-”buttons.

If you have chosen „**CONT**” or „**BURST**” you can set the frequency (**Hz**) and pulse width (**μS**).

If you have chosen „**MOD**” you can set the following parameters:

- upper limit of frequency (**Hz HI**)
- low limit of frequency (**Hz LO**)
- upper limit of pulse width (**μS HI**)
- low limit of pulse width (**μS LO**)
- low limit of intensity (**AdJ**) - the display in the right corner below shows values from 0.5 to 1.0. This means that the intensity of current can increase and decrease within the modulation time between 50% (=0.5) and 100% (=1.0) of the originally set mA-value. The upper limit of the intensity (100%) is always the mA value adjusted by the user originally. The intensity in a modulated program, can never increase above the 100% of the originally set intensity.
- modulation time (**TM** - possible choice: 2 sec. up to 60 sec)

Programming of Dense-Disperse-Stimulation according to Prof. Han (**CONT + BURST**) is not possible.

If you have chosen "**WORK/REST**"- rhythm you can set the following parameters:

- frequency (**Hz**)
- pulse width (**µS**)
- work time (**WORK**)
- rest time (**REST** = recovery time)
- increase / decrease ramp up time of pulse (**RU**)
- increase / decrease ramp down time of pulse (**Rd**)
- synchronous (**SY**) or alternating (**AL**) output of the channels A and B. If you have selected "**SY**" you can program **delay-time (DL)** from 0 sec. up to 5 sec. between the channels A and B.

6. By pressing the „S”button, you can set the next phase of the program. The active phase is indicated by a flashing numeral at the top of the screen.

The total number of numerals indicates the number of configured phases in this program.

You can create up to 5 phases of your own in the configurable program. If you want to reduce a configured program of 5 phases to less phases, you simply set the time of the first unwanted phase to zero.

7. The programming is completed by pressing the „P” button and the program is stored.

Selectable frequency and pulse width settings:

Continuous-Mode: 1 to 150 Hz, 50 to 350 µS

Burst-Mode: 35 to 150 Hz, 50 to 200 µS

Modulated-Mode: 1 to 150 Hz, 50 to 350 µS

Work/Rest-Mode: 1 to 150 Hz, 50 to 450 µS

If frequency is > 100 Hz, pulse width is limited to 350 µS.

The settings can be made as follows.

- Frequency: 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 14, 15 Hz, and then increasing in 5 Hz steps until max. 150 Hz
- Pulse width: 50 to 450 µS in steps of 10 microseconds
- Intensity: 50-100% in steps of 10%
- Modulation time: 2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 15, 20, 30, 40, 50, 60 seconds
- Work time: 2-99 seconds in steps of 1 second
- Ramp up time: 0,1-9,9 seconds in steps of 0,1 second
- Ramp down time: 0,1-9,9 seconds in steps of 0,1 second
- Rest time: 2-99 seconds in steps of 1 second
- Delay: 0,0-5,0 seconds in steps of 0,1 second
- Session time: 1-599 minutes in steps of 1 minute

Programs

PRG 01 = Muscle rehabilitation, muscle strengthening, increasing resistance to muscular fatigue
- high exertion

PRG 02 = Muscle rehabilitation, muscle strengthening - less exertion

PRG 03 = Relaxation and detoxification of the muscle, acceleration of the removal of toxic
metabolic waste products, prevention of muscle cramps

PRG 04 = Reduction and/or elimination of pain

PRG 05 can be configured in up to 5 different phases


Program	01	02	03	04	05
Mode	W/R	W/R	MOD	HAN	selectable
Frequency	50 Hz	35 Hz	8-12 Hz	100/2 Hz	1-150 Hz
Pulse width	300 μ S	450 μ S	300-200 μ S	200/350 μ S	50-450 μ S
Running Time	20 min.	20 min.	30 min.	30 min.	1-599 min.
Ramp up	1 sec.	1 sec.			0,1-9,9 sec.
Ramp down	1 sec.	1 sec.			0,1-9,9 sec.
Work Period	10 sec.	3 sec.			2-99 sec.
Rest Period	20 sec.	15 sec.			2-99 sec.

Abbreviations:

AdJ:	Adjustment of intensity
AL:	Alternate function of the channels A and B
BURST:	The burst mode delivers 2 blocks of 9 impulses each per second
CONT:	In the continuous mode a permanent unchanged current is used
dA:	Day
dEL:	Delete
dL:	Delay between channels A and B
HI:	Upper value (high)
Hz:	Hertz is the measurement unit for the frequency of impulses per second.
LO:	Lower value
mA:	Milli-amp is the measurement unit for the electrical current intensity.
MOD:	In the modulated mode the frequency and/or the pulse width and/or the intensity changes in a predefined rhythm automatically.
Rd:	Ramp down
RU:	Ramp up
SEC:	Seconds
SY:	Synchronous function of the channels A and B
TI:	Time
TM:	Time of modulation-interval
W/R:	Work/Rest mode
µS:	Microsecond is the measurement unit for the length of every single impulse.


LOCK-button and Statistics

The **Rehalito** offers the ability to lock one program in place so that other programs or settings cannot be selected or changed. This allows the doctor or therapist to ensure that the patient uses only the program that has been prescribed. The therapist can read information about the compliance of the patient on the display at any time he wants. This way it is possible to evaluate, whether the device was used according to the instructions of the doctor (quantity, intensity, etc.) For this purpose the LOCK-button in the battery-compartment has to be pressed with a small pin once (the metal-pins of the electrode-lead wires are well suited). The lock symbol appears on the display and the program cannot be changed by pressing the „P“ -button anymore. Pressing the LOCK-button a second time allows reading the following statistics for the treatment in the display:

If the program was not locked, it is necessary to press the LOCK button twice to read the memory. The statistics are always recorded regardless of whether the unit is locked or unlocked. The statistics should be reset after 60 days of use or if the device is used by another person.



To get more detailed information press the + button of the left channel (A). You will see the data of the first session of the first day.

LOCK-button and Statistics

The next sessions for the selected day are displayed - one after the other - by pressing the „S”- button. The next day can be obtained by pressing the + button of channel A again. Hold A+ or A- for three seconds to return to the global statistics screen.

If you require information about time and date of one single session please press any of the right +/--buttons (channel B). With these buttons you can scroll between session data, date and time.

The data stored in memory is retained after viewing, so that therapy-control is possible anytime. A maximum of 60 days of data can be stored with 5 sessions per day. When 60 days of data are reached, consecutive days will overwrite the earliest days. The day number will continue increasing (61, 62, ... up to max. 180) and will be displayed when the unit is unlocked until the **Rehalito** is reset.

If you press the P button for three seconds, you return to the home screen of the selected program. The device will not be locked. Press the lock button once to lock the device if it is required.

If you do not press „P” for 3 seconds, you can turn off the unit to leave it in lock mode when it is restarted.

The statistics should be reset before a new patient uses the device. Press LOCK to enter the statistics display and press „S” twice for 5 seconds each time. All settings in program PC05 will be reset as well. The **Rehalito** is set to its default state and all statistics are deleted.

Trouble Shooting Guide

Fault	Possible Cause	Solution
Trip off - Electrode symbol flashing, current won't increase above 8 mA)	Electrode lost adhesion to skin.	Check electrodes and replace if necessary.
	Electrode to dry and loses adhesion to skin.	Replace electrode.
	Lead wire not connected properly.	Check lead wire connections and ensure they are fully inserted.
	Lead wire is faulty.	Replace lead wire.
	Device is damaged.	Return device to manufacturer for repair.
Device does not power on.	Flat battery or no batteries installed in device.	Check / Replace batteries.
	Battery installed in wrong orientation.	Remove batteries and install according to diagrams in battery compartment.
	Battery connector damaged or flattened.	Remove batteries, check metal parts in battery compartment for damage.
	Device is damaged.	Return to manufacturer for repair.
Device powers off during use.	Batteries almost flat.	Check battery icon during use and replace batteries if low.
	Device is damaged.	Return to manufacturer for repair.
Current output is unstable	Device is in use near a device emitting strong radio or magnetic interference.	Move at least four meters away from the offending device and resume treatment.
Treatment is painful or uncomfortable when using skin electrodes.	Patients skin has cuts or abrasions in the area of electrode placement.	Discontinue use until skin is fully healed.
Skin becomes red after treatment.	Device is being used for periods which are too long for the patient.	Use device for a shorter treatment duration.
Treatment is uncomfortable, tingling sensations during use.	Intensity too high for required treatment.	Reduce intensity setting until treatment is comfortable. Note: Tingling sensations are normal and should be felt.

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To achieve the best possible effect, one of the electrodes should be put on the motor point of the respective muscle if possible. Try to find out this point by slightly moving the electrode around. The other electrode is fastened on the remote muscle end. However, you also obtain a good effect if you put the electrodes respectively on the opposite muscle ends.

■ = red cable connection and ■ = black cable connection

While stimulating, please keep this device in reach, to be able to finish the treatment at once, if current becomes uncomfortable.

Please ensure you follow the warnings and electrode placement advices on pages 23-25 of this manual.

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In most cases it is very easy to find the suitable points for the placing of the electrodes. It is recommended to seek advice on electrode placement from your doctor or physiotherapist.

We suggest placing the electrodes on either side of the pain if no advice is available.

The distance between the electrodes should not be more than 20 cm. The polarity of the electrodes has no significance in TENS therapy.


Examples of possible placement-points are on the following pages.

The intensity should be set to a level that can be felt and then reduced slightly. The intensity setting should not cause visible muscle contraction.

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
Please ensure you follow the warnings and electrode placement advice on pages 23-25 of this manual.

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
















In accordance with European regulations [1], used electrical and electronic equipment may not be discarded in household or unsorted municipal waste. Use separate collection facilities. Contact your local government or council for information regarding the collection systems available. If electrical appliances are disposed of in landfills or dumps, hazardous substances can leak into the ground water and get into the food chain damaging the health and well being of people and animals.

[1] European Community directive 2012/19/EU on waste electrical and electronic equipment.

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The batteries should be replaced when they are drained or stop functioning. The batteries supplied with the **Rehalito** are recyclable. Remove the batteries from the unit and follow your local recycling guidelines.

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A list of literature regarding the topics of this manual can be found on our website „www.mtrplus.com“.

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
Reserved for changes.

Printing or copying of this manual - also in an abridged version - is forbidden and can only be allowed with a written approval of “MTR+ Vertriebs GmbH”.

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Electromagnetic Compatibility

If the device is used according to the instructions given in this manual, it corresponds to the safety requirements for medical electrical systems IEC 60601-1. **Rehalito** was checked and is in accordance with the limits set for medical products in the standard IEC 60601-1-2:2020. These limits serve to the safety of the user regarding electromagnetic disturbances.

The **Rehalito** is a Class B device and is suitable for home healthcare environments as well as typical clinical environments.

Guidance and manufacturer's declaration - electromagnetic emissions		
The Rehalito unit is intended for use in the electromagnetic environment specified below. The customer or the user of the Rehalito should ensure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Rehalito uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Rehalito is suitable for use in all establishments, including domestic and medical establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not Applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not Applicable	

Please adhere to these warnings when using the device:

Don't use **Rehalito** near active HF surgical equipment or RF shielded rooms used for magnetic resonance imaging, where the intensity of EM disturbances is high.

Use of this equipment adjacent to, or placement on top of other equipment should be avoided in case the output of the **Rehalito** is effected. If this is necessary, Both the **Rehalito** and the other equipment should be monitored to verify that they are both operating normally.

Use of accessories and lead wires other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used not closer than 30 cm (12 inches) to any part of the **Rehalito**, including lead wires specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Electromagnetic Compatibility

Guidance and manufacturer's declaration - electromagnetic Immunity			
The Rehalito is intended for use in the electromagnetic environment specified below. The customer or the user of the Rehalito should assure that it is used in such an environment			
Immunity Test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %
Electrical fast transient/burst IEC 61000-4-4	Power supply lines: ±2 kV input/output lines: ±1 kV	Not Applicable	Rehalito is powered by batteries
Surge IEC 61000-4-5	line(s) to line(s): ±1 kV. line(s) to earth: ±2 kV. 100 kHz repetition frequency	Not Applicable	Rehalito is powered by batteries
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° phases 0% U_T for 0.5 cycle And at single phase 0° 0% U_T for 1 cycle 70% U_T for 25/30 cycles 0% U_T for 250 cycles 0% U_T for 300 cycles	Not Applicable	Rehalito is powered by batteries
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m 50 Hz / 60 Hz	30 A/m 50 Hz / 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Conducted RF IEC61000-4-6	150 KHz to 80 MHz: 3 Vrms 6 Vrms (in ISM and amateur radio bands) 80 % AM at 1 kHz	150 KHz to 80 MHz: 3 Vrms 6 Vrms (in ISM and amateur radio bands) 80 % AM at 1 kHz	Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC61000-4-3	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	
NOTE 1: U_T is the a.c. mains voltage prior to application of the test level. NOTE 2: At 80 MHz the higher frequency range applies. NOTE 3: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

Electromagnetic Compatibility

Guidance and manufacturer's declaration - electromagnetic Immunity

The **Rehalito** is intended for use in the electromagnetic environment specified below. The customer or the user of the **Rehalito** should assure that it is used in such an environment

	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Modulation (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	385	380 – 390	TETRA 400	Pulse modulation 18 Hz	1,8	0,3	27
	450	380 – 390	GMRS 460, FRS 460	FM c) ± 5 kHz deviation 1 kHz sine	2	0,3	28
	710	704 – 787	LTE Band 13, 17	Pulse modulation 217 Hz	0,2	0,3	9
	745						
	780						
	810	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0,3	28
	870						
	930						
	1720	1 700 – 1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0,3	28
	1845						
	1970						
	2450	2 400 – 2 570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0,3	28
	5240	5 100 – 5 800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0,2	0,3	9
	5240						
	5785						

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